ENGLISH

General
3M™ Filtek™ P60 Restorative material is a visible-light activated, radiopaque, restorative composite. It is designed for use in posterior restorations. The filler in 3M Filtek P60 restorative is zirconia/silica. The inorganic filler loading is 61% by volume (without silane treatment) with a particle size range of 0.01 to 3.5 µm. 3M Filtek P60 restorative contains BIS-GMA, UDMA and BIS-EMA resins. A 3M dental adhesive is used to permanently bond the restoration to the tooth structure. The restorative is available in a variety of shades. It is packaged in traditional syringes.

Indications

3M Filtek P60 restorative is indicated for use in:

- Direct posterior restorations
- Core Build-ups
- Spinting
- Indirect restorations including inlays, onlays and veneers

Precautions

3M Filtek P60 restorative contains methacrylates. A small percentage of the population is known to have an allergic response to acrylate resins. To reduce the risk of allergic response, minimize exposure to these materials. In particular, exposure to uncured resin should be avoided. **Use of protective gloves and a no-touch technique is recommended.** If restorative material contacts skin, wash immediately with soap and water. Acrylates may penetrate commonly used gloves. If restorative contacts glove, remove and discard glove, wash hands immediately with soap and water and then reglove. If accidental contact with eyes or prolonged contact with oral soft tissues occurs, flush immediately with large amounts of water.

Instructions for Use

I. Preliminary

A. Prophy: Teeth should be cleaned with pumice and water to remove surface stains.

B. Shade Selection: Before isolating the tooth, select the appropriate shade(s) of restorative material.

C. Isolation: A rubber dam is the preferred method of isolation. Cotton rolls plus an evacuator can also be used.

II. Posterior Restorations

A. Cavity Preparation: Prepare the cavity. Line and point angles should be rounded. No residual amalgam or other base material should be left in the internal form of the preparation that would interfere with light transmission and therefore, the hardening of the restorative material.

B. Pulp Protection: If a pulp exposure has occurred and if the situation warrants a direct pulp capping procedure, place a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ Vitrebond™ Light Cure Glass Ionomer Liner/Base. 3M Vitrebond liner/base may also be used to base areas of deep cavity excavation. See 3M Vitrebond liner/base instructions for details.

C. Placement of Matrix: Place a thin dead-soft metal, or a precontoured mylar or a precontoured metal matrix band and insert wedges firmly. Burnish the matrix band to establish proximal contour and contact area. Adapt the band to seal the gingival area to avoid overhangs.

D. Adhesive System: Follow the manufacturer’s instructions regarding etching, priming, adhesive application and curing.

E. Dispensing the Composite: Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of the restorative when dispensing is completed, turn the handle counterclockwise a half turn to stop paste flow. Immediately replace syringe cap. If not used immediately, the dispensed material should be protected from light.

F. Placement:

1. Using a nonmetallic placement instrument, place restorative into the cavity in increments no thicker than 2.5mm.

Placement hints:

a) To aid in adaptation, the first 1mm layer may be placed and adapted to the proximal box.

b) Avoid intense light in the working field.

c) A condensing instrument (or similar device) can be used to adapt the material to all of the internal cavity aspects.

2. Light cure each increment 20 seconds by exposing its entire surface to a high intensity visible light source, such as a 3M™ Curing Light. Hold the light guide tip as close to the restorative as possible during light exposure.

3. *Slightly* overfill the cavity to permit extension of composite beyond cavity margins. Contour and shape with appropriate composite instruments.

G. Finishing: Contour restoration surfaces with fine finishing diamonds, burs or stones. Contour proximal surfaces with 3M™ Finishing Strips.

H. Adjust Occlusion: Check occlusion with thin articulating paper. Centric and lateral excursion contacts should be examined. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

I. Polishing: Polish with 3M™ Sof-Lex™ Discs and Strips or with white stones and rubber points where discs are not suitable.

III. 3M Filtek P60 restorative Indirect Procedure For Inlays, Onlays Or Veneers

A. Dental Operatory Procedure

1. Shade selection: Choose the appropriate shade(s) of 3M Filtek P60 restorative prior to isolation.

2. Preparation: Prepare the tooth.

3. Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer’s instructions of the impressioning material chosen. Any 3M impressioning system may be used.

B. Laboratory Procedure.

1. Pour the impression of the preparation with die stone. Place pins at the preparation site at this time if a “triple tray” type of impression was used.

2. Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base the cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model to an adequate articulator.

3. If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working cast.

4. Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is being used.

5. Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, then add another thin layer.

6. Add the first third of composite to the floor of the preparation, stay short of the margins, light cure for 20 seconds.

7. Add second third of composite. Allow for the last third (incisal) to include the contact areas, light cure for 20 seconds.

8. Place the die back into the articulated arch, add the last third of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured incisal increment. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process.

9. With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy.

10. Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.

11. Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, then polish.

C. Dental Operatory Procedure

1. Roughen the interior surfaces of the indirect restoration.

2. Clean the prosthesis in a soap solution in an ultrasonic bath. Rinse thoroughly.

3. Cementation: Cement the prosthesis using a 3M resin cement system by following manufacturer’s instructions.

IV. Storage and Use:

A. Do not expose restorative materials to elevated temperatures or intense light.

B. Unopened kits should be kept at room temperature (15°-27° or 59°-80°F).

C. Do not store materials in proximity to eugenol-containing products.

D. The composite pastes are designed for use at room temperature of approximately 21°-24°C or 70°-75°F. Shelf life at room temperature is 3 years.

V. Warranty

3M will replace product that is proven to be defective. 3M does not accept liability for any loss or damage, direct or consequential, arising out of the use or the inability to use these products. Before using, the user should determine the suitability of the product for its intended use and user assumes all risk and liability whatsoever in connection therewith.

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ENGLISH

General
3M™ Filtek™ P60 Restorative material is a visible-light activated, radiopaque, restorative composite. It is designed for use in Posterior restorations. The filler in 3M Filtek P60 restorative is zirconia/silica. The inorganic filler loading is 61% by volume (without silane treatment) with a particle size range of 0.01 to 3.5 µm. 3M Filtek P60 restorative contains BIS-GMA, UDMA and BIS-EMA resins. A 3M dental adhesive is used to permanently bond the restoration to the tooth structure. The restorative is available in a variety of shades. It is packaged in traditional syringes.

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2. Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base the cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model to an adequate accuracy.

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