GENERAL INFORMATION

Scotchbond™ dual cure dental adhesive is a two-part bonding adhesive which is cured by exposure to light from a dental visible light curing unit or will auto cure when covered by 3M restorative and the restorative cured. The adhesive bonds to both dentin and etched enamel and chemically combines with 3M restoratives to provide retention and sealing. It is recommended for all classes of restorations.

PRECAUTIONS FOR DENTAL PERSONNEL AND PATIENTS

Scotchbond dual cure adhesive precautions:

1. Scotchbond dual cure adhesive contains BIS-GMA and TEGDMA resins. A small percentage of the population is known to have an allergic response to acrylate resins. To reduce the risk of allergic response, minimize exposure to these materials. In particular, exposure to uncured resins should be avoided. If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If skin contact occurs, wash skin with soap and water.

ETCHANT PRECAUTIONS

Etchant contains 35% phosphoric acid. Protective eyewear for patients and dental staff is recommended when using etchants. Avoid contact with oral soft tissue, eyes and skin. If accidental contact occurs, flush immediately with large amounts of water. For eye contact, also consult a physician.

ADDITIONAL NOTE

Protect adhesive mixture in the mixing tray from ambient light by appropriately covering the tray.

KEY PLACEMENT STEPS

Prophy.

Prepare conservative cavity.

Bevel Cavosurface margins.

Place pulp protection.

Acid etch cavosurface enamel margins for 15 seconds.

Rinse for 15 seconds and dry thoroughly.

Mix adhesive resin and liquid components.

Apply a single adhesive coating to cavity covering dentin and etched enamel.

Direct air gently across surface to thin adhesive and evaporate solvent.

Light Cure 10 seconds.

Or apply a second adhesive coating and air thin.

Place 3M Restorative.

Cure.

Finish and polish.

Refer to manufacturer's instructions for placement and finishing of restorative materials. Note that if not cured by light exposure after application and before restorative placement, the Scotchbond dual cure dental adhesive will auto cure when covered by the restorative and the restorative cured. If used in this manner, delay contouring and finishing for a minimum of 5 minutes after the restorative has been cured.

ENAMEL ETCHING

On primary teeth or teeth high in fluoride, a 90-120 second etch is advised.

SCOTCHBOND DUAL CURE DENTAL ADHESIVE

1. Do not use Resin or Liquid to adjust viscosity of restorative pastes.

2. Use one well of two-well mixing tray for mixing resin and liquid and one for dispensing the etching gel.

3. Clean resin/liquid mixing tray by first wiping with tissue and then with alcohol or acetone; clean etching gel dispensing tray with water and dry.

4. Use blue brush to apply etchant; use black brush to apply Scotchbond dual cure dental adhesive. Clean the etchant brush with water and dry; clean the adhesive brush with alcohol or acetone immediately after use.

STORAGE AND USE

1. Do not expose materials to elevated temperatures or intense light.

2. Do not store materials in proximity to eugenol-containing products.

3. This system is designed to be used at room temperature approximately 21° - 24° C or 70° - 75° F.

4. Shelf life at room temperature is 36 months. See outer package for expiry date.

WARRANTY

3M will replace such product that is proved to be defective. 3M does not accept liability for any loss or damage, direct or consequential, arising out of the use of or the inability to use these products. Before using,

the user should determine the suitability of the product for its intended use and user assumes all risk and liability whatsoever in connection therewith.

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