## ENGLISH

### DESCRIPTION

Clinpro™ Sealant, manufactured by 3M ESPE, is a light-cure, fluoride releasing pit and fissure sealant with a unique color change feature. Clinpro Sealant is pink when applied to the tooth surface, and changes to an opaque light yellow color when exposed to light. The pink color aids the dental professional in the accuracy and amount of material placed during the sealant procedure. When light cured, the pink sealant will transform to an opaque light yellow color.

Note: The change of color from pink to light yellow is not a cure indicator. Sealant needs to be cured with a dental curing light for the recommended exposure time. A sealant exhibiting any pink coloring is not completely cured.

### TECHNICAL INFORMATION

- Meets ISO 6874 (Dental resin based pit and fissure sealant)
- Meets ANSI/ADA Spec 39 (Pit and fissure sealant)
- BIS-GMA / TEGDMA resin composition
- Unfilled
- Curing light must have minimum output of 400 mW/cm<sup>2</sup>
- Use at room temperature

### PRECAUTIONS FOR DENTAL PERSONNEL AND PATIENTS:

- **Etchant Precautions:** Scotchbond™ Etchant, manufactured by 3M ESPE, contains 35% by weight phosphoric acid. Protective eyewear for patients and dental staff is recommended when using etchants. Avoid contact with oral soft tissue, eyes, and skin. If accidental contact occurs, flush immediately with large amounts of water. For eye contact, immediately rinse with plenty of water and seek medical attention
- **Sealant Precautions:** Clinpro Sealant contains acrylate resins. Avoid use of this product on patients with known acrylate allergies. To reduce the risk of allergic response, minimize exposure to these materials. In particular, exposure to uncured resin should be avoided. **Use of protective gloves and a no-touch technique is recommended.** If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If sealant contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If irritation persists, consult a physician.

### INDICATIONS:



Clinpro Sealant is designed for sealing the enamel pits and fissures of teeth to aid in the prevention of caries.

### DISPENSING SEALANT

Follow the directions corresponding to the dispensing system chosen. Sealant is light sensitive. Exposure to overhead operator lights will initiate the color change and curing.

#### Syringe

1. Protective eyewear is recommended for patients and staff when using a syringe type dispenser.
2. Prepare delivery system: Remove cap from syringe and SAVE. Twist a disposable tip securely onto the syringe. Holding the tip away from the patient and any dental staff express a small amount of material onto a mix pad or 2x2 gauze to assure the delivery system is not clogged. If clogged, remove the tip and express a small amount of material from the syringe. Remove any visible plug, if present, from the syringe opening. Replace syringe-tip and again check flow from tip. If clog remains, discard dispensing tip and replace with a new one.
3. At the completion of the procedure remove used syringe-tip and discard. Twist on storage cap. Storage of the syringe with a used dispensing tip, or without the storage cap will allow drying or curing of the product and consequent clogging of the system. Replace storage cap with a new dispensing tip at next use.
4. Disinfection: Discard used syringe tip and replace with syringe storage cap. Disinfect the capped syringe in the same manner as recommended by the ADA and CDC for non-immersible dental items. Council on Dental Materials, Instruments, and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988).

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
#### Sellador Instrucciones

#### Bottle

1. Dispense 1 to 2 drops of sealant into the mix well. **Immediately** slide cover over well to protect from light.
2. Re-cap sealant bottle.
3. After removing material from well always replace cover slide.
4. Disinfection: Disinfect the bottle following procedures for non-immersible dental items as stated under "Syringe #4".

Disinfect mix well and applicator handles following disinfecting solution manufacturer's recommendations.

#### STORAGE AND USE:

- Replace caps on syringes and bottles immediately after use.
- Do not expose materials to elevated temperature.
- Do not store materials in proximity to eugenol-containing products.
- The etchant and sealant is designed for use at room temperatures of approximately 21-24°C or 70-75°F
- Shelf life at room temperature is 24 months. See expiration date on sealant package.

#### APPLICATION GUIDE:

The acid etch technique requires care, particularly for isolation and prevention of contamination. The enamel to be bonded must be cleaned, and thoroughly washed and

dried, and maintained free from contamination after etching procedure and prior to sealant placement.

As an alternative to the acid etch technique the use of Adper™ Prompt™ Self Etch Adhesive, manufactured by 3M ESPE, may be substituted. Follow the Adper Prompt Adhesive instructions, section titled: Instructions for bonding light-cured dental sealants. If Adper Prompt Adhesive is used, enamel must still be thoroughly cleaned and isolated, follow instructions 1-4 in the acid etch technique.

### ACID ETCH TECHNIQUE:

1. **Check air/water syringe.** Blow a jet of air from syringe onto a glove or mirror. If small droplets are seen the syringe must be adjusted so only air is expressed. Any moisture contamination during certain stages of this procedure will compromise the integrity of a sealant.
2. **Select Teeth.** Teeth must be sufficiently erupted so that a dry field can be maintained. The morphology of the pits and fissures should be deep.
3. **Clean Enamel.** Thoroughly clean teeth to remove plaque and debris from enamel surfaces and fissures. Rinse thoroughly with water.

**Note:** Do not use any cleaning medium that may contain oils. If cleaning teeth using an air-polisher that utilizes sodium bicarbonate, it is recommended to repeat steps 5 and 6, or apply hydrogen peroxide to the surface for 10 seconds to neutralize the sodium bicarbonate, then rinse thoroughly with water.

4. **Isolate Teeth and Dry.** While a rubber dam provides the best isolation, cotton rolls used in conjunction with isolation shields, are acceptable. Use saliva ejection device, or high volume evacuation if possible.
5. **Etch Enamel.** Apply a generous amount of etchant to all enamel surfaces to be sealed, extending beyond the anticipated margin of the sealant. Etch for a minimum of 15 seconds, but no longer than 60 seconds.
6. **Rinse Etched Enamel.** Thoroughly rinse teeth with air/water spray to remove etchant. Remove rinse water with suction. Do not allow patient to swallow or rinse. If saliva contacts the etched surfaces, re-etch for 5 seconds and rinse.
7. **Dry Etched Enamel.** Thoroughly dry the etched surfaces. Air should be oil and water free. The dry etched surfaces should appear as a matte frosty white. If not, repeat steps 5 and 6. **DO NOT ALLOW THE ETCHED SURFACE TO BE CONTAMINATED.** Clinical studies have clearly shown that moisture contamination of these surfaces is the main cause for failure of pit and fissure sealants. Immediately apply sealant.
8. **Apply Sealant.** Using the syringe needle tip, or a brush slowly introduce sealant into the pits and fissures. Do not let sealant flow beyond the etched surfaces. Stirring the sealant with the syringe-tip during or after placement will help eliminate any possible bubbles, and enhance the flow into the pit and fissures. An explorer may also be used. Cure the sealant by exposing it to light from a 3M ESPE light curing unit, or other curing unit of comparable intensity. A 20-second exposure is needed for each surface. The tip of the light should be held as closely as possible to the sealant, without actually touching the sealant. When set, the sealant forms a hard, opaque film light yellow in color with a slight surface inhibition.

**9. Evaluate Sealant.** Inspect sealant for complete coverage and voids. If surface has not been contaminated, additional sealant may be added. If contamination has occurred re-etch, rinse, and dry prior to placing more sealant.

10. **Dismissal.** Wipe the sealant with a cotton applicator to remove the thin sticky film on the surface. Check occlusion and adjust as required.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

#### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

#### Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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3M, ESPE, Clinpro, Scotchbond, Adper and Prompt are trademarks of 3M or 3M ESPE AG. ©3M 2004. All rights reserved.