Precautions for Dental Personnel and Patients:**Etchant Precaution:**

Etching Gel contains 35 weight % phosphoric acid. Etching Liquid contains 37 weight % phosphoric acid.

Protective eyewear for patients and dental staff is recommended when using etchants. Avoid contact with oral soft tissue, eyes and skin. If accidental contact occurs, flush immediately with large amounts of water. For eye contact, also consult a physician.

Instructions for Use:

Scotchbond etching gel and 3M etching liquid are intended for use in the acid etch technique. Certain products require an etch of enamel only, while other products require etching enamel and dentin. Please etch according to the instructions of the Scotchbond adhesive that is used.

1. Clean Enamel: If enamel to be etched is not freshly cut, clean with a pumice-water slurry. Rinse well. Prophyl pastes containing oil or fluoride should not be used as these additives may interfere with etching.

2. Isolate: A dry operating field is essential for the acid etch procedure. Rubber dam is the preferred method of isolation.

3. Etching: Dispense the etching gel or etching liquid. Using an appropriate applicator instrument, apply the gel or liquid to enamel for **15 seconds**.

4. Rinse: Thoroughly rinse with oil-free water for 15 seconds. Remove rinse water with suction. Do not allow patient to rinse. If saliva contacts the etched surfaces, re-etch for 5 seconds and rinse. Re-isolate the etched surfaces if using cotton rolls.

5. Dry: Thoroughly dry the etched surfaces. Air should be oil and water free. The dry etched enamel should appear frosty white. If not, repeat steps two through five.

DO NOT ALLOW THE ETCHED SURFACES TO BE CONTAMINATED BY ANYTHING. If contamination occurs, re-etch, rinse and dry as above.

Immediately proceed with the bonding procedure.

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Storage and Use:

- Etching gel and etching liquid are designed for use at room temperatures of approximately 70 - 75° F (21 - 24° C).
- Shelf life at room temperature is 36 months. See outer package for expiry date.

Warranty

3M will replace such product that is proved to be defective. 3M does not accept liability for any loss or damage, direct or consequential, arising out of the use of or the inability to use these products. Before using, the user should determine the suitability of the product for its intended use and user assumes all risk and liability whatsoever in connection therewith.

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