


3M Animal Care Equine Solutions



CLINTON ANDERSON
DOWNUNDER HORSEMANSHIP®
Sponsor



Animalintex® Poultice Pad or Hoof Pad



Ready-to-use poultice — just add water

Animalintex Poultice is a convenient, ready-to-use poultice that comes in a pad, not a bucket. This unique poultice contains a mild antiseptic (boric acid) and Basson (Tragacanth), a natural poulticing agent to promote faster healing.

Key Features:

- Helps draw out infected wounds and reduce inflammation
- Aids in wound healing and relieves minor stiffness and soreness
- Absorbs exudates; provides moist, clean healing site; cushions sore areas
- Retains moisture/warmth and helps prevent external contamination due to its occlusive backing
- Available in two sizes: pre-cut hoof shape or 8" x 16" pad (can be cut to any size).

Recommended Uses:

- Hot poultice — for abscesses, boils, infected and dirty wounds, cracked heels, thorns, laminitis, seedy toe, corns
- Cold poultice — strains, sprains, sore shins, splints, bruises, capped hocks and elbows
- Can also apply to hoof after training to help cool down/condition hoof
- Dry dressing — over open wounds, pressure padding

"Animalintex is a very effective poultice for foot abscesses and sole bruises. It works well for strains and sprains of the legs, too. Because the poultice is contained within the bandage material, it is easy to use with virtually no mess."

Dr. Steven White,
Sport Horse Veterinary Service, Gretna, NE

Gamgee® Highly Absorbent Padding



More than just a cotton roll

Gamgee Highly Absorbent Padding is a hospital quality 100% cotton roll with a non-woven cover that is used for a wound dressing, protective padding, support and warmth/insulation.

Key Features:

- Non-woven, low adherent cover helps prevent bandage from sticking to wound and helps keep cotton fibers out of the wound site
- Minimizes wound disruption and provides a fiber and particle barrier
- Highly absorbent
- Strong (even when wet) protective padding
- Autoclavable: can be cut to size and sterilized (steam or EO)
- Two sizes:
 - 12" wide x 11.5' (knee/hock to fetlock or foal applications)
 - 18" wide x 7.5' (knee/hock to pastern or draft horse applications)

Recommended Uses:

- Surgical/wound dressing over open wounds and for infected and discharging wounds
- Limb immobilization
- Pressure padding to help reduce or prevent edema. Supports and protects fetlocks, hocks, tendons, absorbs sweat, helps prevent pressure sores
- Warmth/insulation — helps keep tendons and muscles warm when exercising, eventing, traveling
- Stall and shipping wraps — economical, durable, can be re-used
- First aid kit

"The Gamgee cotton wrap is easy to handle, conforms to the leg well and provides excellent uniformity and density to a bandage."

Dr. Kevin Voller
Anoka Equine Veterinary Services, Elk River, MN

Activate Activated Carbon Dressing



Helps promote wound healing

Activated carbon dressing promotes rapid wound closure and cleansing by adsorbing exudate and bacteria from the wound into the dressing.

Key Features:

- Multi-layered dressing contains activated carbon in a unique dressing composition
- Low-adherent cover on dressing: Helps prevent sticking to wound
- Helps reduce wound odor
- Helps establish clean granulation bed for healing
- Individual, sterile dressing that can be left on for up to 3 days

Recommended Uses:

- Heavily infected and chronic wounds — works to adsorb bacteria from wound into dressing
- Can be used as a first aid dressing
- Useful as an aid to manage difficult to heal wounds

3M™ Rundown Patches



Protection for performance horse's fetlock point

3M Rundown Patches are made from tough, flexible, polyurethane film.

Key Features:

- One-sided adhesive circle that is applied directly to 3M Vetrap Bandaging Tape or 3M Equisport Equine Support Bandage (never directly to the skin)

Recommended Uses:

- Protection for horse's fetlock point from rundown burns

3M™ Vetrap™ Bandaging Tape



The professional's choice in self-adherent wraps

Layers "bond" together to provide support bandages that stay put, allowing better fit and lower risk of bandage slipping.

Key Features:

- Bandage sticks to itself (not to hair)
- Conforms snugly around contours that are difficult to bandage
- Versatile. High-performance bandaging tape for a broad spectrum of applications – support, compression and bandage secural

Recommended Uses:

- Supports sprains, splints and convalescent injuries
- Protects wounds and holds sterile dressings in place
- Ideal overwrap for Gamgee padding or Animalintex poultice
- Tail wrap

3M™ Equisport™ Equine Support Bandage



Unparalleled support and protection for performance horses

Self-adherent elastic wrap for high support. Intended as a racing support bandage for professional trainers.

Key Features:

- Absorbs up to 90% more impact force¹ than Vetrap tape without compromising range of motion
- 4.5 times the elastic content of 3M Vetrap Bandaging Tape
- University tested, race track proven
- Energy-absorbing bandage

Recommended Uses:

- Safeguards legs from hypertension damage
- Provides protection from severe rundown
- Useful for horses returning from injury
- Hoof wrap to secure Animalintex Hoof Poultice

1. Equisport bandage splint figure 8 bandage compared to traditional Vetrap tape figure 8 track bandage



Clinton Anderson, host of "Downunder Horsemanship," is a leading trainer, clinician and the first person to win the prestigious Road to the Horse Colt Starting Competition twice in a row.

On 3M™ Vetrap™ Bandaging Tape:

"There are a handful of items no barn should be without. For me, 3M Vetrap Bandaging Tape is something I just have to have at hand at all times. Whether I'm treating injuries, or preventing them, I've found that no self-adherent wrap works better. You'll never find a better product, mate."

On Animalintex® Poultice Pad:

"The Animalintex Poultice Pad has quickly become one of my favorite timesavers in the barn. Because it's ready to use in a self contained pad — there's no muss, no fuss — it just goes straight to work on my horse's wounds and injuries."

On Gamgee® Highly Absorbent Padding:

"I make sure we have Gamgee Padding on hand at all times, on the road and back at the ranch. It's so versatile — one day you might need it to cover a wound that you don't want anything sticking to and the next day you might need it for a stall or shipping wrap. You won't want to be without it either, mate!"

For product questions, call the 3M Customer Service Helpline at 1-800-848-0829.



**Medical Division
Animal Care Products**

3M Center, Building 275-5W-05
St. Paul, MN 55144-1000
U.S.A.
1 800 848-0829
www.3M.com
www.3M.com/animalcare

Animalintex and Gamgee are trademarks of Robinson Healthcare Ltd.
3M, Vetrap and Equisport are trademarks of 3M.
Refer to product packages for complete application instructions.

Please recycle. Printed in U.S.A.
© 3M 2008. All rights reserved.
70-2009-2691-6